


Protocol for new spacer cloning in pCas9



Vector Digest

-Digest 1-2ug of pCas9 with *Bsa*I (NEB)

x ul	pCas9
1 ul	<i>Bsa</i> I (NEB)
5 ul	10X NEB Buffer
0.5 ul	100X BSA
y ul	ddH ₂ O
<hr/>	
50 ul	

-Gel purify digested pCas9 (important for successful cloning)

Oligo Design

-In order to clone a spacer sequence into pCas9, synthesize two oligos (IDT) of the following form:

5' - AACNNNNNNNNNNNNNNNNNNNNNNNNNNNNNG - 3'
 ||||| ||||| ||||| ||||| ||||| |||||
 3' -NNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNCAAA - 5'

Phosphorylation

1 ul	oligo I (100 uM)
1 ul	oligo II (100 uM)
5 ul	10X T4 Ligase buffer (NEB)
1 ul	T4 PNK (NEB)
42 ul	ddH ₂ O
<hr/>	
50 ul	

Annealing

-Add 2.5 ul of 1M NaCl to the phosphorylated oligo pairs.

-Incubate 5' @ 95 and slowly cool down to room temperature (use a thermocycler).

(Alternatively, use a heat block and take the block out of the heater and let cool naturally for 2 hours)

-Dilute annealed oligos 10 times

Ligation

1 ul	gel purified, BsaI digested pCas9
1 ul	diluted annealed oligos
2 ul	10X T4 Ligase buffer (NEB)
1 ul	T4 ligase
15 ul	ddH ₂ O
<hr/>	
20 ul	

-Incubate at RT for 2h or 16C for O/N

Transformation

-*E.coli* DH5 α