



Protocol: PCR of sgRNAs for Illumina sequencing

MATERIALS

- Ex Taq DNA polymerase (Clontech, RR001A)
- PCR plates
- P5 & P7 primers (listed at the end)
- 70% EtOH
- AMPure purification system (Beckman Coulter, 63880)
- 96-well round bottom plate (Costar, 07-200-103)
- Magnet (e.g. Alpaqua, A0011322)

PROTOCOL

PCR set-up

Preferably, prepare mix inside a PCR hood, after cleaning the surface with DNase Away and 70% EtOH. If possible, we recommend setting up at least 4 parallel PCRs for a given sample.

Final contents of each reaction:

- 10 μ L 10x reaction buffer
- 8 μ L dNTP
- 0.5 μ L P5 primer mix, 100 μ M
- 1.5 μ L ExTaq polymerase
- up to 10 μ g of genomic DNA or 200 ng of plasmid DNA
- 10 μ L of P7 primer 5 μ M
- up to 100 μ L with water

1. Make a master mix of reaction buffer, dNTP, P5 primer mix, taq polymerase, and water. Aliquot into a PCR plate.
2. Add template DNA.
3. Add a unique P7 primer to barcode each individual reaction.

Thermal cycler parameters

Wait for block to reach 95°C before adding samples.

1. 95°C, 1 minute
2. 95°C 30 seconds (denaturation)
3. 53°C 30 seconds (annealing)
4. 72°C 30 seconds (extension)
Back to step 2, total of 28 cycles
5. 72°C 10 minutes
6. 4°C forever

AMPure XP-PCR purification (recommended)

1. Pool PCR products into an eppendorf (15-30 μ L per well is typically sufficient).
2. Distribute 100 μ L of pooled products to a 96-well round bottom plate
3. Resuspend the magnetic beads included in the AMPure XP reagent, add 100 μ L of beads to each well
4. Mix thoroughly by pipetting up and down 5 times, try not to make bubbles, and incubate at room temperature for 5 minutes. This step binds PCR products 100bp and larger to the magnetic beads. The color of the mixture should appear homogenous after mixing.
5. Place the reaction plate onto a magnet for 2 minutes to separate beads from the solution. Wait for the solution to clear or you see a brown ring around the perimeter of the well before proceeding to the next step.
6. Aspirate the cleared solution from the reaction plate and discard. This step must be performed while the reaction plate is situated on the magnet. Do not disturb the ring of separated magnetic beads.
7. Add 200 μ L of 70% ethanol to each well and incubate for 30 seconds at room temperature; aspirate the ethanol and discard.
8. Repeat step 7 once more for a total of two ethanol washes.
9. Remove the plate from the magnet and dry plate for 1 minute and no longer than 4 minutes. A longer dry time (the bead ring appears cracked) will significantly decrease elution efficiency.
10. Add 50 μ L of TE buffer to elute the PCR product (elution is rapid—approximately 30 seconds).
11. Place the plate back onto the magnet for ~ 2 minutes.
12. Remove the eluted product and store in an eppendorf. The sample is now ready to be sequenced.

Alternatively, standard gel extraction can be performed. If so, we recommend to isopropanol precipitate the sample post-extraction to remove contaminating guanidine isothiocyanate that is commonly used in gel extraction protocols, which obscures quantitation due to its absorbance at 230 nm.

PCR PRIMERS

A mix of P5 primers with stagger regions of different length is necessary to maintain sequence diversity across the flow-cell. A minimum of 8 primers is recommended. Typically, we mix all 8 primers together in one master mix of P5 primers and use that same mix in all wells.

P5/P7 flowcell attachment sequence

Illumina sequencing primer

Vector primer binding sequence

Stagger region / Barcode region

P5 primer for either vector:

5'**AATGATACGGCGACCACCGAGATCTACACTCTTTCCCTACACGACGCTCTTCCGATCT**[s]**TGTGGAAAGGACGAAACACCG**

P7 primer for use with lentiGuide (product size = 354 nt)

5'**CAAGCAGAAGACGGCATAACGAGATNNNNNNNNGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCTTCTACTATTCTTTCCCTGCACTGT**

P7 primer for use with lentiCRISPRv2 (product size = 285 nt)

5'**CAAGCAGAAGACGGCATAACGAGATNNNNNNNNGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCTCCAATTCCCACTCCTTTCAAGACCT**

| P5 primers | |
|-----------------|--|
| Name | Sequence |
| P5 0 nt stagger | AATGATACGGCGACCACCGAGATCTACACTCTTTCCCTACACGACGCTCTTCCGATCTTTGTGGAAAGGACGAAACACCG |
| P5 1nt stagger | AATGATACGGCGACCACCGAGATCTACACTCTTTCCCTACACGACGCTCTTCCGATCTCTTGTGGAAAGGACGAAACACCG |
| P5 2 nt stagger | AATGATACGGCGACCACCGAGATCTACACTCTTTCCCTACACGACGCTCTTCCGATCTGCTTGTGGAAAGGACGAAACACCG |
| P5 3 nt stagger | AATGATACGGCGACCACCGAGATCTACACTCTTTCCCTACACGACGCTCTTCCGATCTAGCTTGTGGAAAGGACGAAACACCG |
| P5 4 nt stagger | AATGATACGGCGACCACCGAGATCTACACTCTTTCCCTACACGACGCTCTTCCGATCTCAACTTGTGGAAAGGACGAAACACCG |
| P5 6 nt stagger | AATGATACGGCGACCACCGAGATCTACACTCTTTCCCTACACGACGCTCTTCCGATCTTGCACCTTGTGGAAAGGACGAAACACCG |
| P5 7nt stagger | AATGATACGGCGACCACCGAGATCTACACTCTTTCCCTACACGACGCTCTTCCGATCTACGCAACTTGTGGAAAGGACGAAACACCG |
| P5 8nt stagger | AATGATACGGCGACCACCGAGATCTACACTCTTTCCCTACACGACGCTCTTCCGATCTGAAGACCCTTGTGGAAAGGACGAAACACCG |

P7 barcode (index) sequences for all designs (NNNNNNNN)

| Well | Sequence to include in P7 primer, 5'-3' | Index read (rev comp of sequence) | Well | Sequence to include in P7 primer, 5'-3' | Index read (rev comp of sequence) |
|------|---|-----------------------------------|------|---|-----------------------------------|
| A01 | CGGTTCAA | TTGAACCG | E01 | TAACTCAA | TTGAGTTA |
| A02 | GCTGGATT | AATCCAGC | E02 | CGTGAGCC | GGCTCACG |
| A03 | TAACTCGG | CCGAGTTA | E03 | ATCAGAGG | CCTCTGAT |
| A04 | TAACAGTT | AACTGTTA | E04 | TATGGAGG | CCTCCATA |
| A05 | ATACTCAA | TTGAGTAT | E05 | GCGTTCAA | TTGAACGC |
| A06 | GCTGAGAA | TTCTCAGC | E06 | CGCAAGAA | TTCTTGCG |
| A07 | ATTGGAGG | CCTCCAAT | E07 | CGACAGCC | GGCTGTCC |
| A08 | TAGTCTAA | TTAGACTA | E08 | CGACTCGG | CCGAGTCG |
| A09 | CGGTGACC | GGTCACCG | E09 | TACAAGAA | TTCTTGTA |
| A10 | TACAGAGG | CCTCTGTA | E10 | CGCAGATT | AATCTGCG |
| A11 | ATTGTCAA | TTGACAAT | E11 | ATTGCTCC | GGAGCAAT |
| A12 | TATGTCTT | AAGACATA | E12 | GCACTCGG | CCGAGTGC |
| B01 | ATTGGATT | AATCCAAT | F01 | ATGTTCTT | AAGAACAT |
| B02 | ATACTCGG | CCGAGTAT | F02 | ATGTCTCC | GGAGACAT |
| B03 | TATGAGAA | TTCTCATA | F03 | GCACTCAA | TTGAGTGC |
| B04 | GCACAGTT | AACTGTGC | F04 | TAGTAGCC | GGCTACTA |
| B05 | CGTGGATT | AATCCACG | F05 | CGTGTCAA | TTGACACG |
| B06 | TAGTAGAA | TTCTACTA | F06 | GCGTTCTT | AAGAACGC |
| B07 | GCACGATT | AATCGTGC | F07 | GCCAAGCC | GGCTTGCC |
| B08 | CGGTAGCC | GGCTACCG | F08 | GCACCTCC | GGAGGTGC |
| B09 | TAGTTCTT | AAGAACTA | F09 | GCACCTGG | CCAGGTGC |
| B10 | TACAAGTT | AACTTGTA | F10 | GCCAGACC | GGTCTGGC |
| B11 | ATCACTGG | CCAGTGAT | F11 | CGCAAGCC | GGCTTGCC |
| B12 | CGCATCAA | TTGATGCG | F12 | TACATCAA | TTGATGTA |

| | | | | | |
|-----|----------|-----------|-----|----------|-----------|
| C01 | GCACGACC | GGTCGTGC | G01 | GCGTAGCC | GGCTACGC |
| C02 | TACACTCC | GGAGTGTA | G02 | CGACAGAA | TTCTGTCCG |
| C03 | CGGTCTAA | TTAGACCG | G03 | TAGTCTGG | CCAGACTA |
| C04 | ATGTTCCG | CCGAACAT | G04 | ATCAAGTT | AACTTGAT |
| C05 | CGTGGACC | GGTCCACG | G05 | TAGTAGTT | AACTACTA |
| C06 | ATTGAGCC | GGCTCAAT | G06 | ATACTCTT | AAGAGTAT |
| C07 | TAGTTCGG | CCGAACTA | G07 | CGGTAGTT | AACTACCG |
| C08 | CGGTGAGG | CCTCACCG | G08 | ATACGAGG | CCTCGTAT |
| C09 | CGTGAGTT | AACTCACG | G09 | CGCACTGG | CCAGTGCG |
| C10 | ATCAGATT | AATCTGAT | G10 | TACAGACC | GGTCTGTA |
| C11 | TAGTGATT | AATCACTA | G11 | GCGTGACC | GGTCACGC |
| C12 | CGGTTCGG | CCGAACCG | G12 | TATGTCGG | CCGACATA |
| D01 | TATGGACC | GGTCCATA | H01 | CGACTCTT | AAGAGTCG |
| D02 | GCCAAGTT | AACTTGGC | H02 | GCGTTCGG | CCGAACGC |
| D03 | CGCAGACC | GGTCTGCG | H03 | ATACCTAA | TTAGGTAT |
| D04 | CGACCTCC | GGAGGTCG | H04 | CGGTGATT | AATCACCG |
| D05 | GCCACTGG | CCAGTGGC | H05 | TAACGACC | GGTCGTTA |
| D06 | GCGTAGTT | AACTACGC | H06 | ATACAGCC | GGCTGTAT |
| D07 | CGCAAGTT | AACTTGCG | H07 | CGACGACC | GGTCGTCCG |
| D08 | CGACAGTT | AACTGTCCG | H08 | ATCACTAA | TTAGTGAT |
| D09 | CGCATCTT | AAGATGCG | H09 | CGACCTGG | CCAGGTCCG |
| D10 | ATGTTCAA | TTGAACAT | H10 | TATGTCAA | TTGACATA |
| D11 | GCGTAGAA | TTCTACGC | H11 | TAACCTAA | TTAGGTTA |
| D12 | ATGTAGCC | GGCTACAT | H12 | GCCATCTT | AAGATGGC |