

Platinum Gate TALEN construction protocol (Yamamoto lab)

Ver. 1.0

Jan. 2014

Tetsushi Sakuma, Ph. D.

E-mail: tetsushi-sakuma@hiroshima-u.ac.jp

Takashi Yamamoto Lab.

Department of Mathematical and Life Sciences

Graduate School of Science

Hiroshima University

1-3-1 Kagamiyama, Higashi-Hiroshima, 739-8526 Hiroshima, Japan

TEL: +81-82-424-7446

E-mail: tybig@hiroshima-u.ac.jp

This is a protocol for the construction of highly-active **Platinum TALEN** described in Sakuma *et al.*, Sci Rep, 2013 (<http://dx.doi.org/10.1038/srep03379>). Since the assembly method is based on the “Golden Gate TALEN and TAL Effector Kit” (Addgene; cat#1000000024) and the “Yamamoto Lab TALEN Accessory Pack” (Addgene; cat#1000000030), NAR paper by Cermak *et al.* (<http://dx.doi.org/10.1093/nar/gkr218>) and Genes to Cells paper by Sakuma *et al.* (<http://dx.doi.org/10.1111/gtc.12037>) may be of good help to understand more about the principle of the Platinum Gate system.

Note that this kit is **NOT compatible** with the Golden Gate TALEN and TAL Effector Kit, and does **NOT contain** the reporter vector for single-strand annealing assay. If you need the vector for the evaluation of TALEN activities, you should separately obtain the pGL4-SSA vector (Addgene plasmid ID: 42962) aside from the Platinum Gate TALEN kit.

■ Table of Contents

Materials	• • • • •	2
Methods		
1. TALEN design	• • • • •	3
2. Platinum Gate assembly -STEP 1-	• • • • •	3
3. Platinum Gate assembly -STEP 2-	• • • • •	4
Appendix	• • • • •	5

■ Materials

○ Plasmids

Platinum Gate TALEN kit (Addgene)

- Module plasmids: **p1-4HD, p1-4NG, p1-4NI, p1-4NN**
(16 plasmids, Ampicillin resistant)
- Intermediate array plasmids: **pFUS2_a1a, a2a, a2b, a3a, a3b, a4a, a4b, b1, b2, b3, b4**
(11 plasmids, Spectinomycin resistant)
- Destination vectors: **ptCMV-136/63-VR-HD, NG, NI, NN**
ptCMV-153/47-VR-HD, NG, NI, NN
(8 plasmids, Ampicillin resistant)

○ Reagents (excluding standard PCR enzymes and reagents for bacterial culture)

Maker	Product name	Cat. No.	Size/Unit
NEB	Quick Ligation Kit	M2200S	30 reactions
		M2200L	150 reactions
NEB	BsaI-HF	R3535S	1000 units
		R3535L	5000 units
Thermo Scientific	Esp3I	ER0452	1000 units
NEB	10×T4 DNA Ligase Reaction Buffer	B0202S	6 ml
Life Technologies	ChargeSwitch-Pro Plasmid Mini Kit	CS30250	250 reactions

○ Primers

pCR8_F1: TTGATGCCTGGCAGTTCCT
pCR8_R1: CGAACCGAACAGGCTTATGT
TALE-F: GAGCACCCCTCAACCTGACCCAG
TALE-R: CTCGAAAGCTGGGCCACGATTG

■ Methods

1. TALEN design

Design TALEN target sequence using “TALEN Targeter”.


TALEN Targeter

<https://tale-nt.cac.cornell.edu/>

According to our past results, it is preferable that spacer length is around 15 bp (12-16 bp) and 17 bp (15-19 bp) for 153/47 and 136/63 scaffolds, respectively. Optimal number of repeat array is 15-20. Set “Upstream Base” as “T only”.

2. Platinum Gate assembly -STEP 1-

Assemble RVD repeats for the chosen target sequences. Using our methods, you can make TALENs bearing 6-21 repeats. It is necessary that we choose appropriate combinations of vectors for the desired module numbers.



DOI: 10.1038/srep03379

Note that a3c, a4c and a4d are just expediential names. Use a2b as a3c/a4d, and a3b as a4c, respectively.

To construct each pFUS2_axx_TALE plasmid harboring four modules, perform Golden Gate assembly (see Appendix for detailed reaction conditions).

Transform 0.5-1 μ l of the reaction product directly to chemical competent *E. coli* such as XL1-Blue and streak the transformant on Spectinomycin/X-gal/IPTG LB plate as described in Appendix. Culture O/N at 37°C.

Select the desired clones using colony PCR. Detailed reaction condition is described in Appendix.

Culture the desired clones with 100 ng/ μ l Spectinomycin/LB media. Purify the plasmids using Miniprep kit.

3. Platinum Gate assembly -STEP 2-

Perform 2nd-step assembly using the 4-module ligated plasmids constructed in the section 2, following a protocol described in Appendix.

Transform the reaction product to XL1-Blue and streak the transformant on Ampicillin/X-gal/IPTG LB plate. Culture O/N at 37°C.

Screen for correctly assembled clones by colony PCR. You can use TALE-F/TALE-R primers both for ptCMV-153/47-VR and ptCMV-136/63-VR vectors. The detailed reaction condition is described in Appendix.

Culture the desired clones and purify the plasmids using Miniprep kit.

Platinum Gate TALEN construction protocol

-Appendix-

Bench protocol for STEP1 & 2 (Yamamoto lab)

STEP1	(μ l)				
■ Golden Gate reaction		<u>1 module</u>	<u>2 modules</u>	<u>3 modules</u>	<u>4 modules</u>
25ng/ μ l pFUS2 vector	0.3				
50ng/ μ l module	0.3×1	0.3×2	0.3×3	0.3×4	
10 \times T4 DNA ligase buf.	0.2				
BsaI-HF	0.1	premix \rightarrow 0.4 μ l each			
Quick ligase	0.1				
SDW	1	0.7	0.4	(0.1)	
Total	2				

37 $^{\circ}$ C	16 $^{\circ}$ C	4 $^{\circ}$ C	
5min.	10min.	∞	\rightarrow Additional digestion:
$\times 3$			add 10\timesNEBuf.4 0.25 μl, 10\timesBSA 0.25 μl, BsaI-HF 0.1 μl.
			premix \rightarrow 0.6 μ l each 50 $^{\circ}$ C, 30min. / 80 $^{\circ}$ C, 5min. / 4 $^{\circ}$ C, ∞

■ Transformation

Product	0.5	
XL1-Blue (homemade)	10	on ice >5min. \rightarrow 42 $^{\circ}$ C, 30 sec. \rightarrow on ice \rightarrow plating (100 μ g/ml Spec plate w/ X-gal+IPTG)

■ Colony PCR

10X buf. for Hybripol	0.8						
dNTPs	0.64						
10 μ M primer	0.16+0.16	(pCR8_F1/R1)					
Hybripol DNA polymerase (BIOLINE)*	0.04						
MgCl ₂	0.24	95 $^{\circ}$ C	95 $^{\circ}$ C	55 $^{\circ}$ C	72 $^{\circ}$ C	72 $^{\circ}$ C	4 $^{\circ}$ C
SDW	5.96	30s.	15s.	15s.	15s.	15s.	∞
Total	8			$\times 27$			

*Any other standard Taq polymerase can be used.

STEP2**■Golden Gate reaction** a1a a2a-a2b a3a-a3c a4a-a4d

50ng/μl pFUS2_a 0.6 0.6×2 0.6×3 0.6×4

50ng/μl pFUS2_b 0.6

50ng/μl ptCMV vector 0.3

10×T4 DNA ligase buf.	0.4
Esp3I	0.2
Quick ligase	0.2

premix → 0.8μl each

37°C	16°C	4°C
5min.	10min.	∞

SDW 1.7 1.1 0.5 0

×6

Total 4

→ Additional digestion: add 10×Tango buf. 0.5 μl, 10mM DTT 0.5 μl, Esp3I 0.2 μl.

premix → 1.2μl each 37°C, 1hr / 80°C, 5min. / 4°C, ∞

■Transformation

Product 2

XL1-Blue (homemade) 20

on ice >5min. → 42°C, 30 sec. → on ice → plating

(100μg/ml Amp plate w/ X-gal+IPTG)

■Colony PCR

10X buf. for Hybripol 0.8

dNTPs 0.64

10μM primer 0.16+0.16 (TALE-F/R)

Hybripol DNA polymerase (BIOLINE)* 0.04

MgCl₂ 0.24 95°C | 95°C 66°C 72°C | 72°C 4°CSDW 5.96 30s. | 15s. 15s. 50s. | 50s. ∞

Total 8 ×27

*Any other standard Taq polymerase can be used.